QUICK START GUIDE





THE FIRST-EVER 3-IN-1 CPR MANIKIN

- Train CPR on adult female, adult male, and child
- 2020 AHA guideline compliant CPR feedback




ASSEMBLY




- 1. Remove the chest skin
- 2. Remove the face skin
- 3. Make sure the head is correctly attached to the two buffers linking it to the torso




- 4. Move the lung bag below the jaw
- 5. Connect the mouth and airway system to the lung bag
- 6. Make sure the mouth and airway system is correctly attached to the head



- 7. Please ensure:
 - a. The lung bag is attached on the marked circle using the sticker
 - b. The lung bag rests loosely on the neck



- 8. Fit the face skin, checking that the nose, mouth, and upper tab are in the correct position
- 9. Place either female or male chest skin and attach to the studs



10. Place head locking clips in place

TROUBLESHOOTING TIP





CPR Taylor®'s fully functional nose requires that the face skin (and specifically the nostril skin) is properly lined up. Failure to properly line up the nostril skin within the nose piece (and the mouth skin with the mouthpiece) will result in improper airway functioning.

If the chest stops rising during a practice session, first ensure the mouth and airway system is properly connected to the lung bag and the head is positioned correctly. Check that the lung bag is not moist. If the lung bag is moist, replace it with a new one.

MONITOR GUIDE



Adult | Child Option

Standard & Blind Modes
View the test report 10 seconds after
the session has finished in both modes