



LIFELIKE INFANT MANIKIN


Realistic size & movements
ideal for training infant CPR

INCLUDES

- 1 infant manikin
- 1 wireless monitor
- 1 user manual
- 4 lung bags
- 2 mouthpieces
- 1 transport bag

*9V batteries included with
all manikins.

3-year warranty



REALISTIC BREATHING

and practice oropharyngeal
airway (OPA) placement

HEAD TILT

CPR TRAINING

Wireless monitor provides
2020 AHA guideline
compliant CPR feedback

CHOKING PRACTICE

AED TRAINING

Pads do not leave residue

REALISTIC MOVEMENT

Fully movable arms & legs


COMPRESSION RATE

POST-SESSION TEST REPORT ON COMPRESSION RATE

CHEST RECOIL


CHEST COMPRESSION

POST-SESSION TEST REPORT ON COMPRESSION DEPTH & CHEST RECOIL



STANDARD & BLIND MODES

View the test report 10 seconds after
the session has finished in both modes



HIGH-QUALITY MOUTHPIECE

Mouthpiece can be fully
disinfected & reused

2 skin tones



FEATURES



- Realistic head tilt/chin lift
- Incredibly lifelike movement — movable arms and legs
- Perform realistic breathing with BVM
- Oropharyngeal airway placement
- Wireless monitor provides CPR feedback on depth, rate, and recoil
- 2020 AHA guideline compliant
- Lightweight at only 3.42 lbs

**3-YEAR
WARRANTY**

*9V batteries included with all manikins.

HIGHLIGHTS



Wireless CPR Feedback Monitor

The wireless monitor easily pairs with Baby Tyler® to offer CPR feedback on rate, depth, and recoil.

- **Standard Mode:** feedback is visible during practice; view test report after session ends
- **Blind Mode:** feedback is not visible during practice; view test report after session ends
- **Test report after training:** available for standard mode and blind mode

High-Quality Mouthpiece

- Mouthpiece can be fully disinfected & reused



SHOP NOW!

ORDERING INFORMATION

Manikins

#11-208	Baby Tyler® - Light Skin	\$155.99
#11-209	Baby Tyler® - Dark Skin	
Includes:		
<ul style="list-style-type: none"> • 1 Baby Tyler® manikin • 4 lungs and 2 mouthpieces • 1 set of Baby Tyler® clothing • 1 wireless CPR feedback monitor • 1 transport bag • User manual 		
#11-210	Baby Tyler® - 4 Pack - Light Skin	\$598.99
#11-211	Baby Tyler® - 4 Pack - Dark Skin	
#11-224	Baby Tyler® - Diversity Pack - 4 Pack	
Includes:		
<ul style="list-style-type: none"> • 4 Baby Tyler® manikins • 16 lungs and 8 mouthpieces • 4 sets of Baby Tyler® clothing • 4 wireless CPR feedback monitors • 1 transport bag • User manual 		
#11-212	4 & 2 Family Pack - Light Skin	\$995.99
#11-213	4 & 2 Family Pack - Dark Skin	
#11-214	4 & 2 Family Pack - Diversity Pack	
Includes:		
<ul style="list-style-type: none"> • 4 CPR Taylor® adult/child manikins with monitors • 4 male chest skins • 4 female chest skins • 20 adult lungs and 8 mouth & airway systems • 2 Baby Tyler® infant manikins with monitors • 8 infant lungs and 4 mouthpieces • 2 sets of Baby Tyler® clothing • 1 transport bag with wheels • Adult/child and infant user manuals 		
#11-215	Starter Pack - Light Skin	\$339.99
#11-216	Starter Pack - Dark Skin	
Includes:		
<ul style="list-style-type: none"> • 1 CPR Taylor® adult/child manikin with monitor • 1 male chest skin • 1 female chest skin • 5 adult lungs and 2 mouth & airway systems • 1 Baby Tyler® infant manikin with monitor • 4 infant lungs and 2 mouthpieces • 1 set of Baby Tyler® clothing • 1 transport bag • Adult/child and infant user manuals 		

Add-Ons & Replacement Parts

#11-220	Baby Tyler® Lung Bags - 10 Pack	\$16.99
#11-225	Baby Tyler® Mouthpiece - Light Skin - 10 Pack	\$9.99
#11-245	Baby Tyler® Mouthpiece - Dark Skin - 10 Pack	\$9.99
#11-226	Baby Tyler® Face Skin - Light Skin	\$9.99
#11-233	Baby Tyler® Face Skin - Dark Skin	\$9.99
#11-227	Baby Tyler® Chest Skin - Light Skin	\$9.99
#11-235	Baby Tyler® Chest Skin - Dark Skin	\$9.99
#11-246	Baby Tyler® Feedback Monitor	\$79.99
#11-248	Baby Tyler® Clothing	\$24.99

***All prices are in USD**