



in collaboration with



**IN STOCK & READY TO SHIP!**

WorldPoint Products® Baby Tyler® helps you realistically train on a variety of simulated patients. The incredibly lifelike Baby Tyler® replicates the size and movements of a real infant, perfect for training infant CPR. Provide 2020 AHA guideline compliant CPR feedback without relying on smartphones — the wireless CPR feedback monitor provides feedback on compression depth, rate, and recoil in addition to providing a post-session test report.

## FEATURES

- Realistic head tilt/chin lift
- Incredibly lifelike movement — movable arms and legs
- Perform realistic breathing with BVM
- Oropharyngeal airway placement
- Wireless monitor provides CPR feedback on depth, rate, and recoil
- 2020 AHA guideline compliant
- Lightweight at only 3.42 lbs

\*9V batteries included with all manikins.

3-YEAR  
WARRANTY



## HIGHLIGHTS

### Wireless CPR Feedback Monitor

The wireless monitor easily pairs with Baby Tyler® to offer CPR feedback on rate, depth, and recoil.

- **Standard Mode:** feedback is visible during practice; view test report after session ends
- **Blind Mode:** feedback is not visible during practice; view test report after session ends
- **Test report after training:** available for standard mode and blind mode

### High-Quality Mouthpiece

Mouthpiece can be fully disinfected & reused.

## ORDERING INFORMATION

#11-208	Baby Tyler - Light Skin	\$155.99
#11-209	Baby Tyler - Dark Skin	
#11-210	Baby Tyler - 4 Pack - Light Skin	\$598.99
#11-211	Baby Tyler - 4 Pack - Dark Skin	
#11-224	Baby Tyler - Diversity Pack - 4 Pack	
#11-212	4 & 2 Family Pack - Light Skin	\$995.99
#11-213	4 & 2 Family Pack - Dark Skin	
#11-214	4 & 2 Family Pack - Diversity Pack	
#11-215	Starter Pack - Light Skin	\$339.99
#11-216	Starter Pack - Dark Skin	

Visit [worldpoint.com](http://worldpoint.com) for add-ons & replacement parts

## INCLUDES

### Individual Manikin

- 1 Baby Tyler manikin
- 4 lungs and 2 mouthpieces
- 1 set of Baby Tyler clothing
- 1 wireless CPR feedback monitor
- 1 transport bag
- User manual

### 4 Pack

- 4 Baby Tyler manikins
- 16 lungs and 8 mouthpieces
- 4 sets of Baby Tyler clothing
- 4 wireless CPR feedback monitors
- 1 transport bag
- User manual

### 4 & 2 Family Pack

- 4 CPR Taylor adult/child manikins with monitors
- 4 male chest skins
- 4 female chest skins
- 20 adult lungs and 8 mouth & airway systems
- 4 CPR Taylor jackets
- 2 Baby Tyler infant manikins with monitors
- 8 infant lungs and 4 mouthpieces
- 2 sets of Baby Tyler clothing
- 1 transport bag with wheels
- Adult/child and infant user manuals

### Starter Pack

- 1 CPR Taylor adult/child manikin with monitor
- 1 male chest skin
- 1 female chest skin
- 5 adult lungs and 2 mouth & airway systems
- 1 CPR Taylor jacket
- 1 Baby Tyler infant manikin with monitor
- 4 infant lungs and 2 mouthpieces
- 1 set of Baby Tyler clothing
- 1 transport bag
- Adult/child and infant user manuals

