

LAERDAL[®] LITTLE ANNE[®] QCPR



The next generation of our Little Anne[®] manikin includes QCPR measurement and feedback technology that takes the guesswork out of instruction—all from the palm of your hand. Little Anne QCPR works with the QCPR Training App that allows you to monitor and coach up to six students at once, and a Learner app that lets students monitor their own performance.

Little Anne QCPR helps instructors improve CPR training quality and efficiency in their courses, while adding a fun and engaging element for students that promotes learning and motivation.

FEATURES

Little Anne® QCPR is still the same durable, lightweight manikin as the original Little Anne, but now includes QCPR technology that provides:

- Real-time feedback on compressions and ventilations
- A summary with an overall score and improvement suggestion
- Details on compression release, depth and rate, ventilation volume, number of compressions/ventilations/cycles



Little Anne QCPR Upgrade Kit

All existing Little Annes can be quickly and easily upgraded to QCPR with the Little Anne QCPR Upgrade Kit



FEEDBACK OPTIONS

Instructor Functionalities

- Free download on App Store
- Connects to 1-6 manikin(s)
- Shows live and summative feedback of CPR performance
- Provides suggestion for improvement
- Results are saved in the app

Learning Functionalities

- Free download on App Store and Google Play
- Connects to 1 manikin
- Shows live and summative feedback of CPR performance
- Provides Suggestion for improvement
- Results are saved in the app

SkillGuide™

- Provides a tethered *plug and play* option
- Can be used standalone or as a supplement to either of the apps

ORDERING INFORMATION

Manikins

- #10-546 Little Anne QCPR - Light Skin
- #10-547 Little Anne QCPR - Dark Skin
- #10-551 Little Anne QCPR Upgrade Kit
Compatible with all existing Little Anne manikins

Accessories and Consumables

- #10-243 Laerdal® SkillGuide™
- #10-308 Little Anne Airways - 24 Pack
- #10-309 Little Anne Airways - 96 Pack
- #20-140 Manikin Face Shields - 6 Rolls of 36
- #40-147 Manikin Wipes - 50 Pack
- #10-374 Manikin Faces - Light Skin - 6 Pack
- #10-084 Manikin Faces - Dark Skin - 6 Pack