



WorldPoint Products® Baby Tyler® helps you realistically train on a variety of simulated patients. The incredibly lifelike Baby Tyler replicates the size and movements of a real infant, perfect for training infant CPR. Provide 2020 AHA guideline compliant CPR feedback without relying on smartphones—the wireless CPR feedback monitor provides feedback on compression depth, rate, and recoil in addition to providing a post-session test report.

FEATURES

- Realistic head tilt/chin lift
- Incredibly lifelike movement—movable arms and legs
- Perform realistic breathing with BVM
- Oropharyngeal airway placement
- Wireless monitor provides CPR feedback on depth, rate, and recoil
- 2020 AHA guideline compliant
- Lightweight at only 3.42 lbs

*9V batteries included with all manikins.



HIGHLIGHTS

Wireless CPR Feedback Monitor

The wireless monitor easily pairs with Baby Tyler to offer CPR feedback on rate, depth, and recoil.

- **Standard Mode:** feedback is visible during practice; view test report after session ends
- **Blind Mode:** feedback is not visible during practice; view test report after session ends
- **Test report after training:** available for standard mode and blind mode

High-Quality Mouthpiece

Mouthpiece can be fully disinfected & reused.

ORDERING INFORMATION

#11-208	Baby Tyler - Light Skin	\$155.99
#11-209	Baby Tyler - Dark Skin	
#11-210	Baby Tyler - 4 Pack - Light Skin	\$598.99
#11-211	Baby Tyler - 4 Pack - Dark Skin	
#11-224	Baby Tyler - Diversity Pack - 4 Pack	
#11-212	4 & 2 Family Pack - Light Skin	\$995.99
#11-213	4 & 2 Family Pack - Dark Skin	
#11-214	4 & 2 Family Pack - Diversity Pack	
#11-215	Starter Pack - Light Skin	\$339.99
#11-216	Starter Pack - Dark Skin	

Visit worldpoint.com for add-ons & replacement parts

INCLUDES

Individual Manikin

- 1 Baby Tyler manikin
- 4 lungs and 2 mouthpieces
- 1 set of Baby Tyler clothing
- 1 wireless CPR feedback monitor
- 1 transport bag
- User manual

4 Pack

- 4 Baby Tyler manikins
- 16 lungs and 8 mouthpieces
- 4 sets of Baby Tyler clothing
- 4 wireless CPR feedback monitors
- 1 transport bag
- User manual

4 & 2 Family Pack

- 4 CPR Taylor adult/child manikins with monitors
- 4 male chest skins
- 4 female chest skins
- 20 adult lungs and 8 mouth & airway systems
- 4 CPR Taylor jackets
- 2 Baby Tyler infant manikins with monitors
- 8 infant lungs and 4 mouthpieces
- 2 sets of Baby Tyler clothing
- 1 transport bag with wheels
- Adult/child and infant user manuals

Starter Pack

- 1 CPR Taylor adult/child manikin with monitor
- 1 male chest skin
- 1 female chest skin
- 5 adult lungs and 2 mouth & airway systems
- 1 CPR Taylor jacket
- 1 Baby Tyler infant manikin with monitor
- 4 infant lungs and 2 mouthpieces
- 1 set of Baby Tyler clothing
- 1 transport bag
- Adult/child and infant user manuals

