

# LIFE/FORM<sup>®</sup> BASIC BUDDY<sup>®</sup> PLUS



Powered by HeartiSense™

## FEATURES

### Basic Buddy® Plus Powered By HeartiSense™


- Dramatic, realistic, affordable, and 2019 AHA Compliant
- Enables accurate, effective training with apps and sensor kit
- Student and instructor apps (included) provide real-time feedback and feature adult capabilities, as well as a realistic interface
- Basic Buddy® manikin simple to assemble and maintain
- Provides an open airway using the head-tilt/chin-lift method with a visible chest rise when ventilated
- Offers an easy-to-simulate obstructed airway
- Xiphoid process provides anatomical reference point for hand placement and compressions
- Training data is saved digitally on the phone or tablet
- Student and instructor apps run on Apple® iPhone® devices. Apple® iPad® version apps will be available February 2019. Android™ version apps will be available in the first half of 2019. (iPad®, Android™, and Apple® tablet/phone not included).

NOTE: Use only AED training devices. When using foam-bodied manikins for AED training, use only pads designed for foam bodies.

## ORDERING INFORMATION

- #10-810 Basic Buddy® Plus
  - #10-811 Basic Buddy® Plus (5 pack)
  - #10-812 Basic Buddy® Plus (10 pack)
  - #10-814 Basic Buddy® Plus Convenience Pack
  - #10-806 Basic Buddy® Plus Upgrade Kit
  - #10-807 Basic Buddy® Plus Upgrade Kit - 5 pack
- 5-year warranty on manikins  
1-year warranty on add-on kit hardware

## BASIC BUDDY® APP



### Student App

Easy, intuitive real-time training and measurement of CPR Feedback for improving students' skills.



### Instructor App

Able to link up to 6 manikins, well organized monitoring screens, and exporting/printing results.