



BRAYDEN CPR TRAINING MANIKINS

The **Brayden CPR Training Manikin** lets your students see the real point of CPR – getting blood flowing to the brain and other organs. In addition to the features of the Brayden CPR Manikin, the **Brayden Pro CPR Training Manikin** also provides healthcare professionals with quantitative, data-driven feedback. The Brayden Pro sends live data to the Brayden Online App via Bluetooth to provide realtime feedback of CPR performance on up to 6 Brayden Pro Manikins that can be reviewed, analyzed, and easily exported.



FEATURES

- Realistic anatomy including head tilt, chin lift, nose pinch, airways and chest rise/fall
- 3 sets of LED lights give instant feedback:
 - Chest compression indicating lights
 - Blood circulation indicating lights
 - CPR quality indicating lights
- Indicators can be configured to the AHA 2020 CPR Guidelines and parameters can be adjusted in the Brayden CPR Pro App
- Brayden Online App*

- Intuitive app and feedback screens provide real-time CPR feedback, including ventilation performance, hands-off time, and hand position

- Bluetooth® link connects up to 6 manikins to app on an Android* tablet, allowing data to be easily exported

- Provides a timeline of assessed performance measured through 8 key CPR metrics for quantitative analysis and debriefing

- Can be set for different responder skill levels, 1 or 2 responder scenarios, ventilation only, or compression only

*The Brayden Online app will work on any current Android version as long as it is Bluetooth-enabled.



MANIKINS

Adults

- Adult Advanced or Pro Manikin
- 1 spare face skin
- 1 spare lung bag
- 10 face shields
- 6 "C" cell batteries
- Carry case with built-in knee pads
- Instructions for use

Infants

- Infant Advanced or Pro Manikin
- 5 double-ventilation filters
- 10 face shields
- 1 Bodysuit
- 4 "AA" alkaline batteries
- Carry case
- Quick gudie & instructions for use

ORDERING INFORMATION

- #10-788 Brayden Advanced CPR Training Manikin#10-774 Brayden Pro CPR Training Manikin#11-131 Brayden Advanced Baby CPR Training Manikin
- #11-130 Brayden Pro Baby CPR Training Manikin

Please note that the Brayden CPR Training Manikin cannot be upgraded to become a Brayden CPR Pro Training Manikin. In order to use the Brayden CPR Pro App with your manikin, you will need to purchase a Brayden CPR Pro Manikin.