



3-IN-1 MANIKIN

Adult male, adult female, and child CPR training for the price of 1 manikin!

INCLUDES

- 1 manikin
- 1 male chest skin
- 1 female chest skin
- 1 wireless monitor
- 5 lung bags
- 2 mouth & airway systems
- 2 head locking clips
- 1 transport bag

*9V batteries included with all manikins.

1-year warranty

LIFE-LIKE BREATHING SYSTEM

Features a functional nose

JAW THRUST MANEUVER

CPR TRAINING

Wireless monitor provides 2020 AHA guideline compliant CPR feedback

LATEX FREE

AED TRAINING

Pads do not leave residue



HIGH-QUALITY AIRWAY SYSTEM

No need to replace the lung between students — saving you time and money!

- High-quality lung can safely be used until moisture is seen
- Mouth & airway system can be fully disinfected & reused
- Pair with a one-way filter CPR mask and/or BVM for hygienic training

2 SKIN TONES




POST-SESSION TEST REPORT ON COMPRESSION RATE

COMPRESSION RATE

CHEST RECOIL

CHEST COMPRESSION

POST-SESSION TEST REPORT ON COMPRESSION DEPTH & CHEST RECOIL



ADULT | CHILD OPTION STANDARD & BLIND MODES

View the test report 10 seconds after the session has finished in both modes